

Document Name: Master Batch Record: Downstream Purification (Affinity Capture and Anion Exchange)			Page 1	
Number: MBR-DSP-201	Version: 1.0	Effective Date: 01APR2026	Status: Effective	Part No: DS-OBX-201

Record Reproduction Certification

This document is an accurate reproduction of the approved Master Manufacturing Record executed for this lot. Handwritten entries are shown in blue; preprinted form text is black.

Product and Lot Information

Product	OBX-201 (recombinant AAV9)	Drug substance lot	DS-26051
Process step	Affinity capture and anion exchange polish	Manufacturing date	08MAY2026
Affinity resin	POROS CaptureSelect AAVX	Affinity column ID	AFF-COL-02
AEX medium	Strong anion exchange monolith (Q)	AEX column ID	AEX-COL-01
Clarified harvest lot	HAR-26051	Target percent full	NLT 80.0 by AEX-HPLC

Referenced SOPs

SOP Number	Title
SOP-3300	Affinity Capture Chromatography (CaptureSelect AAVX)
SOP-3320	Low-pH Elution and Neutralization
SOP-3400	Anion Exchange Full and Empty Separation
SOP-3410	Fraction Collection and Pooling
SOP-1203	Buffer and Solution Preparation
SOP-1318	In-Process UV, Conductivity, and pH Monitoring
SOP-0905	Bioburden and Endotoxin Sampling

1. EQUIPMENT / INSTRUMENTS			
Description	EIN	Calibration Due	Performed By / Date
Chromatography skid	SKID-DSP-01	30SEP2026	SK 08MAY26
Inline UV detector	UV-231	12OCT2026	SK 08MAY26
Inline conductivity monitor	CON-128	05NOV2026	SK 08MAY26
Inline pH meter	PHM-214	20AUG2026	SK 08MAY26
Fraction collector	FRC-04	18OCT2026	SK 08MAY26
Analytical balance	BAL-021	22NOV2026	SK 08MAY26

2. BUFFER PREPARATION

Buffer	Composition	Lot Number	Prepared By / Date	Verified By / Date
Affinity equilibration	20 mM sodium phosphate, 2 mM MgCl ₂ , pH 5.0	BUF-AEQ-0428	SK 08MAY26	PA 08MAY26
Affinity elution	100 mM sodium citrate, 2 mM MgCl ₂ , pH 3.0	BUF-AEL-0428	SK 08MAY26	PA 08MAY26
Neutralization buffer	1 M Tris base, pH 8.5	BUF-NEU-0428	SK 08MAY26	PA 08MAY26
AEX equilibration	20 mM Bis-Tris-propane, pH 9.0	BUF-QEQ-0428	SK 08MAY26	PA 08MAY26
AEX elution (high salt)	20 mM Bis-Tris-propane, 1 M TMAC, pH 9.0	BUF-QEL-0428	SK 08MAY26	PA 08MAY26

3. AFFINITY CAPTURE CHROMATOGRAPHY

Step	Instruction	Parameter	Actual	Performed By / Date	Verified By / Date
3.1	Record clarified harvest load volume.	Load volume (L). Acceptance: record.	9.8	SK 08MAY26	PA 08MAY26
3.2	Record equilibration conductivity.	Conductivity (mS/cm). Acceptance: 8 to 12.	10.2	SK 08MAY26	PA 08MAY26
3.3	Confirm baseline stable before load.	Baseline stable?	☒ X Yes ☐ No	SK 08MAY26	PA 08MAY26
3.4	Record affinity elution pH.	Elution pH. Acceptance: 2.8 to 3.2.	3.0	SK 08MAY26	PA 08MAY26
3.5	Record low-pH hold time before neutralization.	Hold time (min). Acceptance: not more than 15.	8	SK 08MAY26	PA 08MAY26
3.6	Record neutralized pool pH.	pH. Acceptance: 7.0 to 8.0.	7.6	SK 08MAY26	PA 08MAY26
3.7	Record affinity eluate pool volume.	Volume (L). Acceptance: record.	1.9	SK 08MAY26	PA 08MAY26

4. ANION EXCHANGE FULL AND EMPTY SEPARATION

--	--	--	--	--	--

Step	Instruction	Parameter	Actual	Performed By / Date	Verified By / Date
4.1	Record anion exchange load conductivity.	Conductivity (mS/cm). Acceptance: not more than 5.0 to bind.	4.2	SK 08MAY26	PA 08MAY26
4.2	Record anion exchange load pH.	pH. Acceptance: 8.8 to 9.2.	9.0	SK 08MAY26	PA 08MAY26
4.3	Confirm gradient method loaded.	Method AEX-QA-01. Empty peak about 6.3 mS/cm; full peak about 10.1 mS/cm.	☒ X Yes ☐ No	SK 08MAY26	PA 08MAY26
4.4	Record fraction-collection (pooling) start conductivity.	Conductivity (mS/cm). Acceptance: not less than 8.5 to collect the full-capsid peak and exclude the empty shoulder.	9.0	SK 08MAY26	PA 08MAY26
4.5	Record fraction-collection stop conductivity.	Conductivity (mS/cm). Acceptance: record.	13.4	SK 08MAY26	PA 08MAY26
4.6	Record anion exchange pool volume.	Volume (L). Acceptance: record.	1.6	SK 08MAY26	PA 08MAY26
4.7	Record anion exchange load total vector genomes.	Load VG (e13 vg). Acceptance: record.	5.00	SK 08MAY26	PA 08MAY26
4.8	Record anion exchange pool total vector genomes.	Pool VG (e13 vg). Acceptance: record.	3.30	SK 08MAY26	PA 08MAY26

5. POOL AND STEP RECOVERY

Step	Instruction	Parameter	Actual	Performed By / Date	Verified By / Date
5.1	Calculate anion exchange step recovery.	Step recovery (percent). Acceptance: not less than 55 percent. Formula: Step recovery (percent) = (Pool VG / Load	Pool VG (e13 vg): 3.30 Load VG (e13 vg): 5.00 Result (percent): 66.0	SK 08MAY26	PA 08MAY26

		VG) x 100			
5.2	Record percent full capsids confirmed on the analytical worksheet.	Percent full (percent), reference; confirmed on WS-AEX-26051. Acceptance: not less than 80.0.	92.0	SK 08MAY26	PA 08MAY26
5.3	Confirm purified pool bioburden sample collected.	Bioburden sample collected?	☐ X Yes ☐ No	SK 08MAY26	PA 08MAY26

6. PROCESS COMMENTS

Page / Step	Comment	Performed By / Date	Verified By / Date
Step 3.7	Affinity pool volume initially logged as 1.6 L; single-line corrected to 1.9 L against the balance printout. No process impact.	SK 08MAY26	PA 08MAY26

BATCH RECORD APPROVALS

Manufacturing Review	Quality Assurance Review
F. Castellano 05JUN2026	L. Marchetti 05JUN2026

Manufacturing and QA review signatures are applied at batch disposition.