

|                                                                                          |              |                              |                   |                     |
|------------------------------------------------------------------------------------------|--------------|------------------------------|-------------------|---------------------|
| <b>Document Name: Master Batch Record: Production Bioreactor and Triple Transfection</b> |              |                              | Page 1            |                     |
| Number: MBR-USP-201                                                                      | Version: 1.0 | Effective Date:<br>01APR2026 | Status: Effective | Part No: DS-OBX-201 |

### Record Reproduction Certification

This document is an accurate reproduction of the approved Master Manufacturing Record executed for this lot. Handwritten entries are shown in blue; preprinted form text is black.

### Product and Lot Information

|                          |                                            |                       |                        |
|--------------------------|--------------------------------------------|-----------------------|------------------------|
| Product                  | OBX-201 (recombinant AAV9)                 | Drug substance lot    | DS-26052               |
| Process step             | Production bioreactor, triple transfection | Manufacturing date    | 10JUN2026              |
| Cell line                | Suspension HEK293                          | Working cell bank lot | WCB-293-018            |
| Production bioreactor ID | BRX-200-02                                 | Working volume (mL)   | 10,000                 |
| Transfection reagent     | Linear PEI (1 mg/mL)                       | Target harvest        | 72 h post-transfection |

### Referenced SOPs

| SOP Number | Title                                           |
|------------|-------------------------------------------------|
| SOP-3010   | Working Cell Bank Thaw and Seed Train Expansion |
| SOP-3040   | N-1 Seed Bioreactor Operation                   |
| SOP-3110   | Triple Transfection and Production Culture      |
| SOP-3120   | Production Bioreactor In-Process Monitoring     |
| SOP-3150   | Harvest and Lysis                               |
| SOP-1203   | Buffer and Solution Preparation                 |
| SOP-0905   | Bioburden and Endotoxin Sampling                |

| 1. EQUIPMENT / INSTRUMENTS |            |                 |                     |
|----------------------------|------------|-----------------|---------------------|
| Description                | EIN        | Calibration Due | Performed By / Date |
| Production bioreactor skid | BRX-200-02 | 30SEP2026       | TN 10JUN26          |
| Dissolved oxygen probe     | DO-207     | 12OCT2026       | TN 10JUN26          |
| pH probe                   | PHM-209    | 05NOV2026       | TN 10JUN26          |
| Viable cell counter        | VCC-011    | 20AUG2026       | TN 10JUN26          |
| Analytical balance         | BAL-021    | 22NOV2026       | TN 10JUN26          |
| Biosafety cabinet          | BSC-06     | 15DEC2026       | TN 10JUN26          |

## 2. PLASMID AND REAGENT PREPARATION

| Component            | Identity / Composition         | Lot Number     | Prepared By / Date | Verified By / Date |
|----------------------|--------------------------------|----------------|--------------------|--------------------|
| Rep/Cap plasmid      | pOBX-RepCap (AAV9)             | PL-REPCAP-0417 | TN 10JUN26         | PA 10JUN26         |
| Helper plasmid       | pHelper (Ad E2A, E4, VA RNA)   | PL-HELP-0417   | TN 10JUN26         | PA 10JUN26         |
| Transgene plasmid    | pOBX-ITR-transgene             | PL-TRANS-0417  | TN 10JUN26         | PA 10JUN26         |
| Transfection reagent | Linear PEI, 1 mg/mL            | PEI-0392       | TN 10JUN26         | PA 10JUN26         |
| Production medium    | Chemically defined, serum-free | MED-2208       | TN 10JUN26         | PA 10JUN26         |

## 3. SEED AND INOCULATION

| Step | Instruction                                | Parameter                                                  | Actual | Performed By / Date | Verified By / Date |
|------|--------------------------------------------|------------------------------------------------------------|--------|---------------------|--------------------|
| 3.1  | Record viable cell density at inoculation. | Viable cell density (e6 cells/mL). Acceptance: 0.5 to 1.5. | 0.9    | TN 10JUN26          | PA 10JUN26         |
| 3.2  | Record viability at inoculation.           | Viability (percent). Acceptance: not less than 90.         | 95     | TN 10JUN26          | PA 10JUN26         |
| 3.3  | Confirm working volume.                    | Working volume (mL). Acceptance: 10,000 plus or minus 500. | 10,000 | TN 10JUN26          | PA 10JUN26         |

## 4. TRANSFECTION

| Step | Instruction                                       | Parameter                                                                 | Actual                                                        | Performed By / Date | Verified By / Date |
|------|---------------------------------------------------|---------------------------------------------------------------------------|---------------------------------------------------------------|---------------------|--------------------|
| 4.1  | Record viable cell density at transfection.       | Viable cell density (e6 cells/mL). Acceptance: 1.0 to 2.0.                | 1.7                                                           | TN 10JUN26          | PA 10JUN26         |
| 4.2  | Record culture volume at transfection.            | Culture volume (mL). Acceptance: record.                                  | 10,000                                                        | TN 10JUN26          | PA 10JUN26         |
| 4.3  | Calculate total plasmid DNA for the transfection. | Total DNA (mg). Acceptance: 1.0 microgram/mL per SOP-3110. Formula: Total | DNA per mL (microgram/mL): 1.0<br>Culture volume (mL): 10,000 | TN 10JUN26          | PA 10JUN26         |

|     |                                                          |                                                                                                                |                                                                          |            |            |
|-----|----------------------------------------------------------|----------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------|------------|------------|
|     |                                                          | DNA (mg) = DNA per mL (microgram/mL) x Culture volume (mL) / 1000                                              | Result (mg): 10.0                                                        |            |            |
| 4.4 | Calculate PEI mass from the DNA to PEI ratio.            | PEI mass (mg).<br>Acceptance: DNA to PEI 1:2 by weight per SOP-3110.<br>Formula: PEI (mg) = Total DNA (mg) x 2 | Total DNA (mg): 10.0<br>DNA:PEI ratio: 1:2<br>Result (mg): 20.0          | TN 10JUN26 | PA 10JUN26 |
| 4.5 | Split plasmid DNA by molar ratio.                        | Rep/Cap :<br>Helper :<br>Transgene = 1:1:1 (mg each).                                                          | 3.33 / 3.33 / 3.33                                                       | TN 10JUN26 | PA 10JUN26 |
| 4.6 | Record DNA-PEI complexation time before addition.        | Complexation time (min).<br>Acceptance: 15 to 30 min.                                                          | 20                                                                       | TN 10JUN26 | PA 10JUN26 |
| 4.7 | Confirm addition of transfection complex to the reactor. | Complex added to reactor?                                                                                      | <input checked="" type="checkbox"/> X Yes<br><input type="checkbox"/> No | TN 10JUN26 | PA 10JUN26 |

| 5. PRODUCTION CULTURE |                                   |                                                                 |        |                     |                    |
|-----------------------|-----------------------------------|-----------------------------------------------------------------|--------|---------------------|--------------------|
| Step                  | Instruction                       | Parameter                                                       | Actual | Performed By / Date | Verified By / Date |
| 5.1                   | Record temperature setpoint.      | Temperature (C).<br>Acceptance: 37.0 plus or minus 0.5.         | 37.0   | TN 10JUN26          | PA 10JUN26         |
| 5.2                   | Record dissolved oxygen setpoint. | Dissolved oxygen (percent).<br>Acceptance: 40 plus or minus 10. | 40     | TN 10JUN26          | PA 10JUN26         |
| 5.3                   | Record hour 24 viability.         | Viability (percent).<br>Acceptance: record.                     | 96     | TN 10JUN26          | PA 10JUN26         |
| 5.4                   | Record hour 24 culture pH.        | pH.<br>Acceptance: 7.0 to 7.2.                                  | 6.84   | TN 10JUN26          | PA 10JUN26         |
| 5.5                   | Record hour 48 viability.         | Viability (percent).<br>Acceptance: record.                     | 91     | TN 10JUN26          | PA 10JUN26         |
| 5.6                   | Record hour 48 culture pH.        | pH.<br>Acceptance: 7.0 to 7.2.                                  | 7.05   | TN 10JUN26          | PA 10JUN26         |
| 5.7                   | Record hour 72                    | Viability (percent).                                            | 79     | TN 10JUN26          | PA 10JUN26         |

|  |            |                        |  |  |  |
|--|------------|------------------------|--|--|--|
|  | viability. | Acceptance:<br>record. |  |  |  |
|--|------------|------------------------|--|--|--|

| 6. HARVEST |                                        |                                                       |        |                     |                    |
|------------|----------------------------------------|-------------------------------------------------------|--------|---------------------|--------------------|
| Step       | Instruction                            | Parameter                                             | Actual | Performed By / Date | Verified By / Date |
| 6.1        | Record harvest time post-transfection. | Harvest time (h).<br>Acceptance: 48 to 96, target 72. | 80     | TN 14JUN26          | PA 14JUN26         |
| 6.2        | Record harvest viability.              | Viability (percent).<br>Acceptance: not less than 70. | 68     | TN 14JUN26          | PA 14JUN26         |
| 6.3        | Record harvest working volume.         | Volume (mL).<br>Acceptance: record.                   | 9,850  | TN 14JUN26          | PA 14JUN26         |

| 7. IN-PROCESS SAMPLING |                                               |                                              |                                                                          |                     |                    |
|------------------------|-----------------------------------------------|----------------------------------------------|--------------------------------------------------------------------------|---------------------|--------------------|
| Step                   | Instruction                                   | Parameter                                    | Actual                                                                   | Performed By / Date | Verified By / Date |
| 7.1                    | Record in-process vector genome titer (qPCR). | VG titer (e10 vg/mL).<br>Acceptance: record. | 8.9                                                                      | TN 14JUN26          | PA 14JUN26         |
| 7.2                    | Confirm bioburden sample collected.           | Bioburden sample collected?                  | <input checked="" type="checkbox"/> X Yes<br><input type="checkbox"/> No | TN 14JUN26          | PA 14JUN26         |
| 7.3                    | Confirm mycoplasma sample collected.          | Mycoplasma sample collected?                 | <input checked="" type="checkbox"/> X Yes<br><input type="checkbox"/> No | TN 14JUN26          | PA 14JUN26         |

| 8. PROCESS COMMENTS |                                                                                                                                                                                                                                      |                     |                    |
|---------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|--------------------|
| Page / Step         | Comment                                                                                                                                                                                                                              | Performed By / Date | Verified By / Date |
| Step 4.4            | PEI mass initially calculated as 30.0 mg applying a 1:3 ratio; SOP-3110 specifies 1:2. Error identified at verification before complex preparation and recalculated to 20.0 mg. No PEI added at the incorrect mass. Ref DEV-26-0205. | TN 10JUN26          | PA 10JUN26         |
| Step 5.4            | Hour 24 culture pH 6.84, below the 7.0 to 7.2 range. CO2 sparge reduced; pH recovered to 7.08 by hour 30. Ref DEV-26-0206.                                                                                                           | TN 10JUN26          | PA 10JUN26         |
|                     |                                                                                                                                                                                                                                      |                     |                    |

|              |                                                                                                                                                                                                                                  |            |            |
|--------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|------------|
| Step 6.1/6.2 | Harvest performed at 80 h against the 72 h target after a scheduling delay; harvest viability 68 percent, below the not-less-than-70 percent target. Material retained; in-process titer within expected range. Ref DEV-26-0208. | TN 14JUN26 | PA 14JUN26 |
|--------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|------------|

| BATCH RECORD APPROVALS  |                          |
|-------------------------|--------------------------|
| Manufacturing Review    | Quality Assurance Review |
| F. Castellano 16JUL2026 | L. Marchetti 16JUL2026   |

Manufacturing and QA review signatures are applied at batch disposition.