

Analytical Worksheet: Vector Genome Titer by Droplet Digital PCR

Master Batch Record worksheet. Titer is computed from the reportable dilution. Encapsidated genomes are counted after DNase digestion.

1 SAMPLE AND METHOD

Ref	Field	Entry
1.1	Worksheet ID	WS-VG-26052
1.2	Product / lot	OBX-201 / DS-26052
1.3	Method reference	ddPCR (ITR target), USP <1047>, SOP-3220
1.4	Test date	16JUL2026
1.5	Instrument (EIN)	Droplet reader, EIN DDR-006
1.6	Analyst	J. Reyes

2 SYSTEM SUITABILITY

Ref	Check	Criterion	Result	Meets criterion
2.1	No-template control positive droplets	Not more than 10	2	Yes
2.2	Positive control titer	Within qualified range	In range	Yes
2.3	Total accepted droplets (reportable well)	Not less than 10,000	14,200	Yes
2.4	Replicate agreement, percent CV	Not more than 15 percent	3.1	Yes

3 REACTION AND READING

Ref	Parameter	Value
3.1	Reportable dilution factor (D)	6,400,000
3.2	Positive concentration (C, copies per microliter)	625.0
3.3	Reaction volume (R, microliter)	20
3.4	Template volume (V, microliter)	4

4 CALCULATION AND RESULT

Ref	Parameter	Value
4.1	Titer formula	$T = (C \times R \times D \times 1000) / V$
4.2	Vector genome titer (vg/mL)	2.00×10^{13}
4.3	Acceptance criterion	1.6 to 2.4×10^{13} vg/mL

4.4	Meets criterion	Yes

5 ATTACHMENTS

Ref	Attachment
5.1	ddPCR run report attached

6 WORKSHEET APPROVALS

Ref	Performed By / Date	Verified By / Date
6.1	J. Reyes 16JUL2026	M. Delgado 17JUL2026