

Analytical Worksheet: Vector Genome Titer by Droplet Digital PCR

Master Batch Record worksheet. Titer is computed from the reportable dilution. Encapsidated genomes are counted after DNase digestion.

1 SAMPLE AND METHOD

Ref	Field	Entry
1.1	Worksheet ID	
1.2	Product / lot	
1.3	Method reference	ddPCR (ITR target), USP <1047>, SOP-3220
1.4	Test date	
1.5	Instrument (EIN)	
1.6	Analyst	

2 SYSTEM SUITABILITY

Ref	Check	Criterion	Result	Meets criterion
2.1	No-template control positive droplets	Not more than 10		
2.2	Positive control titer	Within qualified range		
2.3	Total accepted droplets (reportable well)	Not less than 10,000		
2.4	Replicate agreement, percent CV	Not more than 15 percent		

3 REACTION AND READING

Ref	Parameter	Value
3.1	Reportable dilution factor (D)	
3.2	Positive concentration (C, copies per microliter)	
3.3	Reaction volume (R, microliter)	20
3.4	Template volume (V, microliter)	4

4 CALCULATION AND RESULT

Ref	Parameter	Value
4.1	Titer formula	$T = (C \times R \times D \times 1000) / V$

4.2	Vector genome titer (vg/mL)	
4.3	Acceptance criterion	1.6 to 2.4 x 10 ¹³ vg/mL
4.4	Meets criterion	

5 ATTACHMENTS

Ref	Attachment
5.1	

6 WORKSHEET APPROVALS

Ref	Performed By / Date	Verified By / Date
6.1		